


Bend the **STATSPLINT™** lengthwise to give it strength.

- 4.25" x 36" rolled splint
- Waterproof
- X Ray Translucent
- Washable and reusable
- Cuts with scissors to any length
- Packaged with 2" x 5' **STATGRIP™** tape

As a multi-purpose emergency splint, the **STATSPLINT™** can be used to immobilize just about any bone in the body, for example, sugar tong, upper arm, ankle stirrup, leg splint etc. A lengthwise "C" curve makes it strong, a reverse "C" curve makes it stronger, and a "T" bend makes it strongest.



STATSPLINT™
by NeeQuip, Inc.
www.statsplint.com
631-689-5514

Authorized Distributor: